Members













SUPERIOR SHOULDER IMMOBILIZER

Product Code: 40-157

Size


CODE	SIZE	
01	Small	
02	Medium	
03	Large	
04	Extra Large	



List of Raw Material:

- Cotton Stockinette
- PU Foam

Product Label



Instruction For Use:

- Used to immobilize the shoulder joint post fracture.
- The shoulder strap is adjustable to allow the correct positioning of the aram after which the restraint band is fitted around the arm and fastened in front with Velcro fastened in front with velcro fasteners.
- Requires normal laundry washing.
- One size fits all.

HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





Exhibitor's at:





Members













Primary Packing



Secondary Packing

Cardboard box made of recycled paper.

Side 1



Side 2



Side 3



Side 4



HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com Web Site: www.hemcmedical.com





Exhibitor's at:





Members













CERTIFICATE OF ANALYSIS

Name of Product: Superior Shoulder Immobilizer Test Report for Finished Products

Name of Product	Superior Shoulder Immobilizer	Specification and model	40-157	
Manufacturer	HOSPITAL EQUIPMENT MANUFACTURING COMPANY (INDIA)	Sample Type	Finished products	
Sampling Location	Ready-for-inspection warehouse of Company	Sampling Date	June - 2024	
Quantity of Sample	100%	Sample Base Number	100 Pcs	
Basis of Test	ISO 9001 : 2015 ISO 13485 : 2016	Lot No. of Production	S-6/700/1590	
Mfg. Date	June-2024	Exp. Date	NA	
Test Conclusion	This lot of product has proved qualified through comprehensive determination after being tested. Date of Reporting: 07.06.2024			
Remarks				



HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





Exhibitor's at:



