

FICHA TÉCNICA

IMAGEN	DESCRIPTION (PRODUCT NAME)	SPECIFICATION AND COMPONENTS	SIZE	MATERIAL COMPONENT	USO	PICTURE PRIMARY AND SECONDARY	PRIMARY AND SECONDARY PACKING
	Cervical Collar REF:FS-238S	The product is composed of front piece, fixed buckle and back piece, among which front piece and back piece shall be molded by EVA foam plastic piece.	S: 13*19*19cm M:13,5*22*20cm L:14*20*20 cm	EVA foam plastic piece.	The product is mainly used for medical personnel to fix patients' cervical vertebrae.	 	the primary container which is a polyethylene bag and the secondary is a paper box

	MANUAL DE USO
Describe how to use the product. You must go step by step how to use.	Paso 1. abra la cinta de velcro y fíjela al cuello Paso 2. arregla el velcro del collar "



FICHA TÉCNICA

IMAGEN	NOMBRE DEL PRODUCTO	ESPECIFICACIONES Y COMPONENTES	TALLA	MATERIALES	USO	IMAGEN DEL EMPAQUE	EMPAQUE PRIMARIO Y SECUNDARIO
	Cervical Collar REF:FS-238S	El producto se compone de una pieza delantera, una hebilla fija y una pieza trasera, entre las cuales la pieza delantera y la trasera serán moldeadas por una pieza de plástico de espuma EVA.	S: 13*19*19cm M:13,5*22*20cm L:14*20*20 cm	pieza de plástico de espuma EVA (etil vinil acetato).	El producto se utiliza principalmente para que el personal médico arregle la vértebra cervical de los pacientes.	 	el contenedor primario que es una bolsa de polietileno y el secundario es una caja de papel

