

- A line of stockings dedicated to prophylaxis of thromboembolism before, during and after surgery.
- Anti-Embolism medical support hose; Unisex ted style for reclined patient recovery.
- Use post-surgery when inactive or bed ridden to help reduce the danger of thromboembolic disease.
- Helps prevent swelling, edema fluid retention, DVT vein clots, pooling of blood, pulmonary embolism.

ANTI-EMBOLISM STOCKING

STOCKING (13-18mmHg)

MAD

Size	. 5	. IVI	<u> </u>	YL.	AAL
LD	36-42cm	36-42cm	36-42cm	36-42cm	36-42cm
сС	30-36cm	33-39cm	36-42cm	39-45cm	42-48cm
сВ	19-23cm	21-25cm	23-27cm	25-29cm	27-31cm

- **LD** (Heel to the length of 3cm below the knee)
 - cC (Calf Circumstance)
 - cB (Ankle Circumstance)

MAG

•	OIZE)	IVI	_	ΛL	AAL
	LG	69-75cm	69-75cm	69-75cm	69-75cm	69-75cm
-	сG	44-52cm	48-56cm	52-60cm	56-64cm	60-68cm
	сС	30-36cm	33-39cm	36-42cm	39-45cm	42-48cm
	сВ	19-23cm	21-25cm	23-27cm	25-29cm	27-31cm

- **▶ LG** (Heel to the length of 5cm below the thigh root)
 - cG (5cm circumference under the thigh root)
 - cC (Calf Circumstance)
 - **cB** (Ankle Circumstance)

item#	Description
MAD-01	AD-KNEE HIGH STOCKINGS
MAG-01	AG-THIGH HIGH STOCKINGS

Indications:

- A valid aid in preventing thromboembolisms linked to surgical interventions.
- Treatments associated with post-traumatic therapies.
- Bruises, sprains and hydrostatic edemas.

Contraindications:

- · Localized dermatitis, erythema, necrosis.
- · Stripping, pulmonary edema, grafts.
- Vascular disorders of cardiac origin.
- Legs of larger circumference

Features:

- Decreasing compression from the ankle to the thigh.
- Extremely soft and elastic fabric, hypoallergenic (latex free).
- Inspection opening on bottom of the foot.



4