



Huaian Wanja Medical Device Co., Ltd.

Mobile Phone: 13515233016 Fax: 86-517-84519182

E-19225B Technical details

Date: 2019-8-2

SR/NO	Image	DESCRIPTION	REREFENCE ITEM NO	commercial and technical name	characteristics	Material composition	USE	Directions of use and points for attention
1		carbon steel surgical blade	10#	surgical blade	<p>A.The Packing Is Done With Seal By The Compound Material Of Polyester And Aluminium Foil With The Nature Of The Good Airtight-Ness And Bacterium Isolation.</p> <p>B.It Is Sterilized By The Co-60 Gamma Radiation And Can Be Used Directly Without Being Re-Sterilized Before Use.</p> <p>C.It Can Prevent The Cross Infections From The Medical Sources And Be Suitable To The Field Rescue.</p>	Carbon Steel And Stainless Steel	This Product Is Used For The Basic Surgical Surgery For Cutting Soft Tissues.	<p>1、Tearing open the package according to the indicated direction by words on the small package.</p> <p>2、Forbidden to use if the small package is damaged.</p> <p>3、Forbidden to use if the appearance is rusty and cracked.</p>
2			11#					
3			12#					
4			15#					
5			22#					
6			23#					
7			24#					