




Technical details of E-20202

DATE:2020-9-8

SR/NO	ITEM NO.	Image	DESCRIPTION 描述	characteristics 特征	Material composition 材质比例	USE 用途
1	AS1019		mini massager	voltage:DC4.5v 3*AAA batteries(excluded), 8 different massage modes:stroke,massage,acupuncture,manipulation,scrapping, cupping,immunotherapy *10 intensity power *15mins timer	ABS 40%+ PCB50% + metal spring 10%	1.alleviating the symptoms of cervical spondylosis , scapulohumeral periarthritis,arthritis and lumbar muscle strain. 2.It can ease peripheral nerve palsy, myalgia 3.It can promote blood circulation 4.It can remove fatigue ,Improve sleep quality 5.Enhance human immunity.Stay away from sub-health 6. A natural and drug free way to help with pain relief.
2	AS-P06		electrode pads	Tens unit pads have universal compatibility with all machines with snap connection. Flexible, comfortable, easy to apply, sticks on lasting comfort and easy to peel off	pvc 20%+conductive gel 80%	electrode pads for tens unit device
3	AS-P01		electrode pads	Tens unit pads have universal compatibility with all machines with snap connection. Flexible, comfortable, easy to apply, sticks on lasting comfort and easy to peel off	no woven cloth 20%+conductive gel 80%	electrode pads for tens unit device