Members














FINGER COT

Product Code: 40-172



### **List of Raw Material:**

- Aluminium Sheet
- PU Foam

### **Product Label**



## **Instruction For Use:**

- Made from aluminium and lined inside with high-density foam for comfort.
- Can be easily fitted or removed from the finger by applying nominal finger pressure to the splint.
- Fits any finger and one size fits all.

**HEMC House:** 

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





Exhibitor's at:





Members













# **Primary Packing**



## Secondary Packing

Cardboard box made of recycled paper.

Side 1







Side 3



Side 4



#### **HEMC House:**

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com Web Site: www.hemcmedical.com





Exhibitor's at:





Members













#### **CERTIFICATE OF ANALYSIS**

Name of Product: Finger Cot

### **Test Report for Finished Products**

Name of Product	Finger Cot	Specification and model	40-172
Manufacturer	HOSPITAL EQUIPMENT MANUFACTURING COMPANY (INDIA)	Sample Type	Finished products
Sampling Location	Ready-for-inspection warehouse of Company	Sampling Date	June - 2024
Quantity of Sample	100%	Sample Base Number	100 Pcs
Basis of Test	ISO 9001 : 2015 ISO 13485 : 2016	Lot No. of Production	S-6/700/1590
Mfg. Date	June-2024	Exp. Date	NA
Test Conclusion	This lot of product has proved quafter being tested.		eporting: 07.06.2024
Remarks			



**HEMC House:** 

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





Exhibitor's at:



