


COMPRESSION STOCKING


CLASS 1 STOCKING (15-21mmHg)



Medical compression products have a well-documented effect for prevention and treatment of venous diseases. The effect of compression therapy is achieved by compressing the limb with a graduated compression, strongest at the ankle and decreasing going up the leg. The compression helps the venous return by aiding the blood flow up the leg against gravity.

The compression is measured in millimetres of mercury(mmHg), an accurate measurement unit for pressure. The compression is created by very elastic yarns knitted in special machines programmed to get the exact pressure at certain positions of the sock. These positions are derived from common quality standards.

Planet can produce compression stockings according to any standard customers wish to apply.



Pressure Rating	cB Pressure Value
	hPa · mmHg

Planet-Y01 20-28 15-21

Planet-Y02 31-43 23-32

▼ Remark : 1mmHg=1.333hPa

AD					
Size	S	M	L	XL	XXL
LD	36-42cm	36-42cm	36-42cm	36-42cm	36-42cm
cC	30-36cm	33-39cm	36-42cm	39-45cm	42-48cm
cB	19-23cm	21-25cm	23-27cm	25-29cm	27-31cm

▼ LD (Heel to the length of 3cm below the knee)

cC (Calf Circumstance)

cB (Ankle Circumstance)

AT/AG					
Size	S	M	L	XL	XXL
LG	69-75cm	69-75cm	69-75cm	69-75cm	69-75cm
cG	44-52cm	48-56cm	52-60cm	56-64cm	60-68cm
cC	30-36cm	33-39cm	36-42cm	39-45cm	42-48cm
cB	19-23cm	21-25cm	23-27cm	25-29cm	27-31cm

▼ LG (Heel to the length of 5cm below the thigh root)

cG (5cm circumference under the thigh root)

cC (Calf Circumstance)

cB (Ankle Circumstance)

Item#	Description
AD-101	AD-KNEE HIGH (OPEN)
AD-102	AD-KNEE HIGH (CLOSE)
AG-101	AG-THIGH HIGH (OPEN)
AG-102	AG-THIGH HIGH (CLOSE)
AT-101	AT-PANTYNOSE (OPEN)
AT-102	AT-PANTYNOSE (CLOSE)



VASTECH MEDICAL compression stockings use high-quality materials. They help improve blood circulation and relieve physical fatigue, edema, and varicose veins with therapeutic results.

COMPRESSION STOCKING

CLASS 2 STOCKING (23-32mmHg)

Pressure Rating	cB Pressure Value hPa	mmHg
-----------------	--------------------------	------

Planet-Y01 20-28 15-21

Planet-Y02 31-43 23-32

Remark : 1mmHg=1.333hPa

Item#	Description
-------	-------------

AD-201 AD-KNEE HIGH (OPEN)


AD-202 AD-KNEE HIGH (CLOSE)

AG-201 AG-THIGH HIGH (OPEN)

AG-202 AG-THIGH HIGH (CLOSE)

AT-201 AT-PANTYNOSE (OPEN)

AT-202 AT-PANTYNOSE (CLOSE)



COMPRESSION STOCKING

CLASS 3 STOCKING (34-46mmHg)

COMPRESSION STOCKING

CLASS 3 STOCKING (34-46mmHg)

Compression stockings use high-quality materials. They help improve blood circulation and relieve physical fatigue, edema, and varicose veins with therapeutic results.

AD

Size · S · M · L · XL · XXL

LD 36-42cm 36-42cm 36-42cm 36-42cm 36-42cm

cC 30-36cm 33-39cm 36-42cm 39-45cm 42-48cm

cB 19-23cm 21-25cm 23-27cm 25-29cm 27-31cm

LD (Heel to the length of 3cm below the knee)

cC (Calf Circumstance)

cB (Ankle Circumstance)

AT/AG

Size · S · M · L · XL · XXL

LG 69-75cm 69-75cm 69-75cm 69-75cm 69-75cm

cG 44-52cm 48-56cm 52-60cm 56-64cm 60-68cm

cC 30-36cm 33-39cm 36-42cm 39-45cm 42-48cm

cB 19-23cm 21-25cm 23-27cm 25-29cm 27-31cm

LG (Heel to the length of 5cm below the thigh root)

cG (5cm circumference under the thigh root)

cC (Calf Circumstance)

cB (Ankle Circumstance)

Pressure Rating	cB Pressure Value hPa	mmHg
-----------------	--------------------------	------

Planet-Y01 20-28 15-21

Planet-Y02 34-43 23-32

Planet-Y03 45-61 34-46

Remark : 1mmHg=1.333hPa

AD

Size · S · M · L · XL · XXL

LD 36-42cm 36-42cm 36-42cm 36-42cm 36-42cm

cC 30-36cm 33-39cm 36-42cm 39-45cm 42-48cm

cB 19-23cm 21-25cm 23-27cm 25-29cm 27-31cm

LD (Heel to the length of 3cm below the knee)

cC (Calf Circumstance)

cB (Ankle Circumstance)

AG

Size · S · M · L · XL · XXL

LG 69-75cm 69-75cm 69-75cm 69-75cm 69-75cm

cG 44-52cm 48-56cm 52-60cm 56-64cm 60-68cm

cC 30-36cm 33-39cm 36-42cm 39-45cm 42-48cm

cB 19-23cm 21-25cm 23-27cm 25-29cm 27-31cm

LG (Heel to the length of 5cm below the thigh root)

cG (5cm circumference under the thigh root)

cC (Calf Circumstance)

cB (Ankle Circumstance)

Item#	Description
-------	-------------

AD-301 AD-KNEE HIGH (OPEN)

AD-302 AD-KNEE HIGH (CLOSE)

AG-301 AG-THIGH HIGH (OPEN)

AG-302 AG-THIGH HIGH (CLOSE)

compression stockings use high-quality materials. They help improve blood circulation and relieve physical fatigue, edema, and varicose veins with therapeutic results.

