











ELASTIC HINGED OPEN PATELLA

Product Code: 40-253

Size


CODE	SIZE	
01	Small	
02	Medium	
03	Large	
04	Extra Large	



List of Raw Material:

- Cotton Stockinette
- PU Foam

Product Label



Instruction For Use:

- Select a size which will be comfortable upon fastening.
- Place the product around affected area.
- Fasten using the velcro take help of a medical practitioner.

HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





















Primary Packing



Secondary Packing

Cardboard box made of recycled paper.

Side 1



Side 2



Side 3



Side 4



HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com Web Site: www.hemcmedical.com





















CERTIFICATE OF ANALYSIS

Name of Product: Elastic Hinged Open Patella Test Report for Finished Products

Name of Product	Elastic Hinged Open Patella	Specification and model	40-253
Manufacturer	HOSPITAL EQUIPMENT MANUFACTURING COMPANY (INDIA)	Sample Type	Finished products
Sampling Location	Ready-for-inspection warehouse of Company	Sampling Date	June - 2024
Quantity of Sample	100%	Sample Base Number	100 Pcs
Basis of Test	ISO 9001 : 2015 ISO 13485 : 2016	Lot No. of Production	S-6/700/1590
Mfg. Date	June-2024	Exp. Date	NA
Test Conclusion	This lot of product has proved qualifier being tested.	-	eporting: 07.06.2024
Remarks			



HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com

Web Site: www.hemcmedical.com





















HEMC House:

D-313, Sector-63, Noida, U. P. 201 301 (INDIA)

Customer care: +91-120-4088 888, Fax: +91-120-4088 777

Email: info@hemcmedical.com Web Site: www.hemcmedical.com







