

IDEAL MEDICAL INDUSTRIES CO., LTD.

LUTIAN BASE, POYANG INDUSTRY ZONE, SHANG RAO, JIANGXI, CHINA

Tel:+86 574 27835985 Fax:+86 574 27835990

E-mail: tank@idealmedical.net www.idealmedical.net

Ideal Sterile Latex Surgical Gloves

Certificate of Analysis

Powder-free

Product Name		Latex Surgical Gloves Powder-free		Product Standard		EN455-1, -2-3	
Lot No.		20241216		Type, Size		6,6.5,7,7.5,8	
MFG Date		12/2024		EXP Date		12/2029	
Lot Quantity		136000 pairs		Sample Quantity		200 pairs	
Sampling		07-02-25~08-02-25		Test Date		08-02-25~15-02-25	
Sampler		Lixia Zhang		Identifier		Lixia Zhang	
Test Item		Index	Test Level		Test Result	Result	
Before Aging	Tensile Strength	≥12 N	S-2 AQL 4.0		13.4~14.6, 0/13	Pass	
Aged	Tensile Strength	≥9 N	S-2 AQL 4.0		12.0~13.2, 0/13	Pass	
6#	Length	≥260 mm	S-2 AQL 4.0		264~268, 0/13	Pass	
	Width	78±5mm	S-2 AQL 4.0		76~82,0/13	Pass	
	Thickness	≥0.10mm	S-2 AQL 4.0		0.29~0.31, 0/13	Pass	
6.5#	Length	≥265 mm	S-2 AQL 4.0		267~272, 0/13	Pass	
	Width	83±5mm	S-2 AQL 4.0		81~85,0/13	Pass	
	Thickness	≥0.10mm	S-2 AQL 4.0		0.29~0.31, 0/13	Pass	
7#	Length	≥270 mm	S-2 AQL 4.0		273~278, 0/13	Pass	
	Width	89±5mm	S-2 AQL 4.0		86~93,0/13	Pass	
	Thickness	≥0.10mm	S-2 AQL 4.0		0.29~0.31, 0/13	Pass	
7.5#	Length	≥270 mm	S-2 AQL 4.0		275~279, 0/13	Pass	
	Width	95±5mm	S-2 AQL 4.0		95~99, 0/13	Pass	
	Thickness	≥0.10mm	S-2 AQL 4.0		0.29~0.31, 0/13	Pass	
8#	Length	≥275 mm	S-2 AQL 4.0		280~282, 0/13	Pass	
	Width	103±6mm	S-2 AQL 4.0		100~107, 0/13	Pass	
	Thickness	≥0.10mm	S-2 AQL 4.0		0.29~0.31, 0/13	Pass	
Leachable Protein		200µg/dm ²	3 pcs		80.0 µg/dm ²	Pass	
Freedom from holes			G-1 AQL1.5		0/200	Pass	
Labeling			S-2 AQL 1.5		0/13	Pass	
sterilize result			N/A		A	Pass	
Biological Test			Non-Irritation to Rabbit's Skin and Non-sensitization to Guinea Pig's Skin.				
Determinant Compliance			This lot conforms to the requirements established by the specifications of EN455-1, EN455-2 and EN455-3.				

Approved By Xiaoling Hong/ QC Dept.

江西云鹤橡塑有限公司
IDEAL MEDICAL INDUSTRIES CO., LTD.
朱轶骐